



Anti-PRR15L (aa 4-94) polyclonal antibody (DPABH-20044)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PRR15L (proline rich 15-like) is a protein-coding gene. Diseases associated with PRR15L include gastric cancer.
Immunogen	antigen sequence, corresponding to amino acids 4-94 (EIGWWKLTL RKKKSTPKVL YEIPDTYAQT EGDAEPPRPD AGGPNSDFNT RLEKIVDKST KGKHKVKSNS GRFKEKKKVR ATLAENPNLF D) of Human PRR15L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PRR15L proline rich 16-like [Homo sapiens]
Official Symbol	PRR15L
Synonyms	PRR15L; proline rich 15-like; ATAD4; proline-rich protein 15-like protein; ATPase family, AAA domain containing 4;
Entrez Gene ID	79170
Protein Refseq	NP_077296.1
UniProt ID	Q9BU68