



# Anti-CCDC144A (aa 913-1009) polyclonal antibody (DPABH-20048)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | CCDC144A (coiled-coil domain containing 144A) is a protein-coding gene. Diseases associated with CCDC144A include smith magenis syndrome. An important paralog of this gene is POTEF.    |
| <b>Immunogen</b>           | antigen sequence, corresponding to amino acids 913-1009 (MESYRCRLAA AVRDCDQSQT ARDLKLDLQFQ TRQEWVRLHD KMKVDMGSLQ AKNEILSEKL SNAESKINSL QIQLHNTRDA LGRESLILER VQRDLSQ) of Human CCDC144A. |
| <b>Isotype</b>             | IgG                                                                                                                                                                                      |
| <b>Source/Host</b>         | Rabbit                                                                                                                                                                                   |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                    |
| <b>Purification</b>        | Immunogen affinity purified                                                                                                                                                              |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                                             |
| <b>Applications</b>        | IHC-P                                                                                                                                                                                    |
| <b>Format</b>              | Liquid                                                                                                                                                                                   |
| <b>Size</b>                | 100 µl                                                                                                                                                                                   |
| <b>Buffer</b>              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                                                                                                                                            |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                                                                                                       |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.                                                                                                    |

## GENE INFORMATION

|                        |                                                                                           |
|------------------------|-------------------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">CCDC144A coiled-coil domain containing 145A [ Homo sapiens ]</a>              |
| <b>Official Symbol</b> | CCDC144A                                                                                  |
| <b>Synonyms</b>        | CCDC144A; coiled-coil domain containing 144A; coiled-coil domain-containing protein 144A; |
| <b>Entrez Gene ID</b>  | <a href="#">9720</a>                                                                      |
| <b>Protein Refseq</b>  | <a href="#">NP_055510.1</a>                                                               |
| <b>UniProt ID</b>      | <a href="#">A2RUR9</a>                                                                    |