



Anti-CCDC144A (aa 859-947) polyclonal antibody (DPABH-20049)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC144A (coiled-coil domain containing 144A) is a protein-coding gene. Diseases associated with CCDC144A include smith magenis syndrome. An important paralog of this gene is POTEF.
Immunogen	antigen sequence, corresponding to amino acids 859-947 (KQQPTSEATS RCHINLDETQ DSKKKLGQIR SEIDLTEAQE TVPSRCLHLD AENEVLQLQQ TLFSMKAIQK QCETLQKNKK QLKQEVVNL) of Human CCDC144A, isoform 3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CCDC144A coiled-coil domain containing 145A [Homo sapiens]
Official Symbol	CCDC144A
Synonyms	CCDC144A; coiled-coil domain containing 144A; coiled-coil domain-containing protein 144A;
Entrez Gene ID	9720
Protein Refseq	NP_055510.1
UniProt ID	A2RUR9