



Anti-C17ORF58 (aa 13-93) polyclonal antibody (DPABH-20051)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C17orf58 (chromosome 17 open reading frame 58) is a protein-coding gene.
Immunogen	antigen sequence: FRVHMLALDS SSCNKPCPEF KPGSRYIVMG HIYHKRRQLP TALLQVLRGR LRPDGLLRS SSSYVKRFNR KREGQIQGAI H, corresponding to amino acids 13-93 of Human C17orf58.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C17ORF58 chromosome 17 open reading frame 59 [Homo sapiens]
Official Symbol	C17ORF58
Synonyms	C17ORF58; chromosome 17 open reading frame 58; UPF0450 protein C17orf58;
Entrez Gene ID	284018
Protein Refseq	NP_858041.2
UniProt ID	Q2M2W7