



Anti-WDR81 (aa 465-589) polyclonal antibody (DPABH-20052)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	WDR81 (WD repeat domain 81) is a protein-coding gene. Diseases associated with WDR81 include cerebellar ataxia, mental retardation, and dysequilibrium syndrome 2 , and dysequilibrium syndrome.
Immunogen	antigen, corresponding to amino acids 465-589 of Human WDR81 (Q562E7).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	WDR81 WD repeat domain 82 [Homo sapiens]
Official Symbol	WDR81
Synonyms	WDR81; WD repeat domain 81; CAMRQ2; PPP1R166; WD repeat-containing protein 81; protein phosphatase 1, regulatory subunit 166;
Entrez Gene ID	124997
Protein Refseq	NP_001157145.1
UniProt ID	Q562E7
Function	transferase activity, transferring phosphorus-containing groups;