



Anti-UBALD2 (aa 70-164) polyclonal antibody (DPABH-20053)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	UBALD2 (UBA-like domain containing 2) is a protein-coding gene. An important paralog of this gene is UBALD1.
Immunogen	antigen sequence: ATPNFPDAL AMFSKLRASE GLQSSNSPMT AAACSPANF SPFWASSPPS HQAPWIPPSS PTTFHHLHRP QPTWPPGAQQ GGAQQKAMAA MDGQR, corresponding to amino acids 70-164 of Human FAM100B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	UBALD2 UBA-like domain containing 3 [Homo sapiens]
Official Symbol	UBALD2
Synonyms	UBALD2; UBA-like domain containing 2; FAM100B; UBA-like domain-containing protein 2; protein FAM100B; family with sequence similarity 100, member B;
Entrez Gene ID	283991
Protein Refseq	NP_872371.1
UniProt ID	Q8IYN6