



Anti-LYRM9 (aa 15-71) polyclonal antibody (DPABH-20054)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LYRM9 (LYR motif containing 9) is a protein-coding gene. An important paralog of this gene is ENSG00000266202.
Immunogen	antigen sequence: QLYRYLLRCC QQLPTKGIQQ HYKHAVRQSF RVHSDDEDNPE RIQQIIKRAI EDADWIM, corresponding to amino acids 15-71 of Human C17orf108.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LYRM9 LYR motif containing 10 [Homo sapiens]
Official Symbol	LYRM9
Synonyms	LYRM9; LYR motif containing 9; HSD24; C17orf108; LYR motif-containing protein 9; UPF0631 protein C17orf108;
Entrez Gene ID	201229
Protein Refseq	NP_001070148.1
UniProt ID	A8MSI8