



Anti-KIAA0753 (aa 850-961) polyclonal antibody (DPABH-20057)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KIAA0753 (KIAA0753) is a protein-coding gene.
Immunogen	antigen sequence, corresponding to amino acids 850-961 of Human KIAA0753. (Q2KHM9)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [KIAA0753 KIAA0754 \[Homo sapiens \]](#)

Official Symbol	KIAA0753
Synonyms	KIAA0753; uncharacterized protein KIAA0753;
Entrez Gene ID	9851
Protein Refseq	NP_055619.2
UniProt ID	Q2KHM9