



Anti-R3HCC1 (aa 55-154) polyclonal antibody (DPABH-20061)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	R3HCC1 (R3H domain and coiled-coil containing 1) is a protein-coding gene. GO annotations related to this gene include nucleotide binding and nucleic acid binding. An important paralog of this gene is R3HCC1L.
Immunogen	antigen sequence: ENDFVHRIQE ELDRFLLQKQ LSKVLLFPPL SSRLRYLIHR TAENFDLLSS FSVGEGWKRR TVICHQDIRV PSSDGLSGPC RAPASCPSRY HGPRPISNQG, corresponding to amino acids 55-154 of Human R3HCC1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	R3HCC1 R3H domain and coiled-coil containing 2 [Homo sapiens]
Official Symbol	R3HCC1
Synonyms	R3HCC1; R3H domain and coiled-coil containing 1; R3H and coiled-coil domain-containing protein 1;
Entrez Gene ID	203069
Protein Refseq	NP_001129580.1
UniProt ID	Q9Y3T6
Function	nucleic acid binding; nucleotide binding;