



Anti-MFSD6L (aa 103-183) polyclonal antibody (DPABH-20063)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	MFSD6L (major facilitator superfamily domain containing 6-like) is a protein-coding gene. An important paralog of this gene is MFSD6.
Immunogen	Recombinant Protein Epitope Signature Tag: NRVHFPCNGS SGLTSTDALP GVTLPVNITS AQESASSHPA KRTAEVEMPG FRNPPGESDR ETFRDLHVYL APSVEGARTT S, corresponding to amino acids 103 - 183 of Human MFSD6L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	MFSD6L major facilitator superfamily domain containing 7-like [Homo sapiens]
Official Symbol	MFSD6L
Synonyms	MFSD6L; major facilitator superfamily domain containing 6-like; major facilitator superfamily domain-containing protein 6-like;
Entrez Gene ID	162387
Protein Refseq	NP_689812.3
UniProt ID	Q8IWD5