



Anti-LRRC45 (aa 294-388) polyclonal antibody (DPABH-20072)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LRRC45 (leucine rich repeat containing 45) is a protein-coding gene.
Immunogen	antigen sequence: RHSIINALKA KLQMTEAALA LSEQKAQDLG ELLATAEQEQ LSLSQRQAKE LKLEQQEAAE RESKLLRDLS AANEKNLLQ NQVDELERKF RCQQE, corresponding to amino acids 294-388 of Human LRRC45
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LRRC45 leucine rich repeat containing 46 [Homo sapiens]
Official Symbol	LRRC45
Synonyms	LRRC45; leucine rich repeat containing 45; leucine-rich repeat-containing protein 45;
Entrez Gene ID	201255
Protein Refseq	NP_659436.1
UniProt ID	Q96CN5