



Anti-CCDC182 (aa 77-153) polyclonal antibody (DPABH-20074)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen sequence: EDAFCEMNGA LVQQEEQAAR VRQLREEED RGIVRNKVLTFLLPREKQLR EHCKRLEDLL LDRGRDALRA TKKSQAD, corresponding to amino acids 77-153 of Human ENSP00000299415
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CCDC182 coiled-coil domain containing 183 [Homo sapiens]
------------------	--

Official Symbol	CCDC182
Synonyms	CCDC182; coiled-coil domain containing 182; coiled-coil domain-containing protein 182; coiled-coil domain-containing protein ENSP00000299415;
Entrez Gene ID	101927581
Protein Refseq	NP_001269473.1
UniProt ID	A6NF36
Function	molecular_function;