



Anti-ZNF347 (aa 47-117) polyclonal antibody (DPABH-20075)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	antigen sequence: LGISCFDLSI ISMLEQGKEP FTLESQVQIA GNPDGWEWIK AVITALSSEF VMKDLLHKGK SNTGEVFQTV M, corresponding to amino acids 47-117 of Human ZNF347
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF347 zinc finger protein 348 \[Homo sapiens \]](#)

Official Symbol	ZNF347
Synonyms	ZNF347; zinc finger protein 347; ZNF1111; zinc finger 1111; zinc finger protein 1111;
Entrez Gene ID	84671
Protein Refseq	NP_001166145.1
UniProt ID	A0A024R4L7
Pathway	Gene Expression.
Function	DNA binding; metal ion binding;