



Anti-EFCAB1 (aa 107-206) polyclonal antibody (DPABH-20078)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen sequence, corresponding to amino acids 107-206 (YCFEVFDLNG DGFISKEEMF HMLKNSLLKQ PSEEDPDEGI) of Human EFCAB1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	EFCAB1 EF-hand calcium binding domain 2 [Homo sapiens]
Official Symbol	EFCAB1

Synonyms	EFCAB1; EF-hand calcium binding domain 1; EF-hand calcium-binding domain-containing protein 1;
Entrez Gene ID	79645
Protein Refseq	NP_001136329.1
UniProt ID	A0A024R7U7
Function	calcium ion binding;