



Anti-LRBA (aa 907-1038) polyclonal antibody (DPABH-20084)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in coupling signal transduction and vesicle trafficking to enable polarized secretion and/or membrane deposition of immune effector molecules.
Immunogen	antigen sequence, corresponding to amino acids 907-1038 of Human LRBA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [LRBA LPS-responsive vesicle trafficking, beach and anchor containing \[Homo sapiens \]](#)

Official Symbol	LRBA
Synonyms	LRBA; LPS-responsive vesicle trafficking, beach and anchor containing; BGL; LBA; CDC4L; CVID8; LAB300; lipopolysaccharide-responsive and beige-like anchor protein; CDC4-like protein; vesicle trafficking, beach and anchor containing;
Entrez Gene ID	987
Protein Refseq	NP_001186211.2
UniProt ID	P50851
Function	molecular_function;