



Anti-ZNF653 (aa 119-203) polyclonal antibody (DPABH-20087)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Transcriptional repressor. May repress NR5A1, PPARG, NR1H3, NR4A2, ESR1 and NR3C1 transcriptional activity.
Immunogen	antigen sequence: NVNCLKNVVI WYEDHKHRCP YEPHLAELDP TFGLYTTAVW QCEAGHRYFQ DLHSPLKPLS DSDPDSDKVG NGLVAGSSDS SSSGS, corresponding to amino acids 119-203 of Human ZNF653.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZNF653 zinc finger protein 654 [Homo sapiens]
Official Symbol	ZNF653
Synonyms	ZNF653; zinc finger protein 653; ZIP67; E430039K05Rik; zinc finger protein Zip67; 67 kDa zinc finger protein;
Entrez Gene ID	115950
Protein Refseq	NP_620138.2
UniProt ID	Q96CK0
Function	DNA binding; metal ion binding;