



Anti-FBF1 (aa 39-115) polyclonal antibody (DPABH-20089)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Keratin-binding protein required for epithelial cell polarization. Involved in apical junction complex (AJC) assembly via its interaction with PARD3. Required for ciliogenesis.
Immunogen	antigen sequence: ARTKSLLGDD VFSTMAGLEE ADAEVSGISE ADPQALLQAM KDLGMDADI LGLKKSNSAP SKKAAKDPGK GELPNHP, corresponding to amino acids 39-115 of Human FBF-1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	FBF1 Fas (TNFRSF6) binding factor 2 [Homo sapiens]
Official Symbol	FBF1
Synonyms	FBF1; Fas (TNFRSF6) binding factor 1; Alb; FBF-1; fas-binding factor 1; albatross; protein albatross; Fas binding protein 1;
Entrez Gene ID	85302
Protein Refseq	NP_001074011.1
UniProt ID	Q8TES7
Function	protein binding;