



Anti-LYNX1 (aa 23-97) polyclonal antibody (DPABH-20091)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Seems to modulate nicotinic acetylcholine receptors.
Immunogen	antigen sequence: IWCHQCTGFG GCSHGSRCLR DSTHCVTTAT RVLSNTEDLP LVTKMCHIGC PDIPSLGLGP YVSIACCQTS LCNHD, corresponding to amino acids 23-97 of Human LYNX1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LYNX1 Ly6/neurotoxin 2 [Homo sapiens]
Official Symbol	LYNX1
Synonyms	LYNX1; Ly6/neurotoxin 1; SLURP2; ly-6/neurotoxin-like protein 1; SLURP-2; Ly-6 neurotoxin-like protein 1; secreted Ly6/uPAR related protein 2; secreted Ly-6/uPAR-related protein 2; secreted Ly-6/uPAR domain-containing protein 2;
Entrez Gene ID	66004
Protein Refseq	NP_076435.1
UniProt ID	Q9BZG9