



Anti-RUNDC1 (aa 354-461) polyclonal antibody (DPABH-20092)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role as p53/TP53 inhibitor and thus may have oncogenic activity.
Immunogen	antigen sequence, corresponding to amino acids 354-461 of Human RUNDC1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [RUNDC1 RUN domain containing 2 \[Homo sapiens \]](#)

Official Symbol	RUNDC1
Synonyms	RUNDC1; RUN domain containing 1; RUN domain-containing protein 1;
Entrez Gene ID	146923
Protein Refseq	NP_775102.2
UniProt ID	Q96C34