



Anti-C1ORF112 (aa 694-774) polyclonal antibody (DPABH-20093)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C1orf112 (chromosome 1 open reading frame 112) is a protein-coding gene. Diseases associated with C1orf112 include essential hypertension , and hypertension.
Immunogen	antigen sequence: SCMFALLAD RSWLLEQHTL EAFTQFAEGT NHEEIVPQCL SSEETKNKVV SFLEKTGFVD ETEAAKVERV KQEKGIWFEP F, corresponding to amino acids 694-774 of Human C1orf112.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C1ORF112 chromosome 1 open reading frame 113 [Homo sapiens]
Official Symbol	C1ORF112
Synonyms	C1ORF112; chromosome 1 open reading frame 112; uncharacterized protein C1orf112;
Entrez Gene ID	55732
Protein Refseq	NP_060656.2
UniProt ID	A0A024R922