



Anti-ZNF48 (aa 27-111) polyclonal antibody (DPABH-20095)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	antigen sequence: GLGSENVISQ PNEFEHTPQE DDLGFKEEDL APDHEVG NVSLKPEGIQNWD DLWVQREGLG KPQPRDRGPR LLGEPRWGQA SSDRA, corresponding to amino acids 27-111 of Human ZNF48 (Q96MX3).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZNF48 zinc finger protein 49 [Homo sapiens]
Official Symbol	ZNF48
Synonyms	ZNF48; zinc finger protein 48; ZNF553; zinc finger protein 553;
Entrez Gene ID	197407
Protein Refseq	NP_001201835.1
UniProt ID	Q96MX3
Function	DNA binding; metal ion binding;