



# Rabbit anti-Human FAM205A polyclonal antibody (DPABH-20105)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	FAM205A (family with sequence similarity 205, member A) is a protein-coding gene.
Immunogen	antigen sequence, corresponding to amino acids 803-921 of Human FAM205A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

Gene Name	<a href="#">FAM205A family with sequence similarity 206, member A [ Homo sapiens ]</a>
-----------	--

<b>Official Symbol</b>	FAM205A
<b>Synonyms</b>	FAM205A; family with sequence similarity 205, member A; C9orf144B; protein FAM205A; transmembrane protein C9orf144B;
<b>Entrez Gene ID</b>	<a href="#">259308</a>
<b>Protein Refseq</b>	<a href="#">NP_001135389.1</a>
<b>UniProt ID</b>	<a href="#">B4DY95</a>