



Anti-CCDC129 (aa 894-1041) polyclonal antibody (DPABH-20108)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC129 (coiled-coil domain containing 129) is a protein-coding gene. An important paralog of this gene is TESPA1.
Immunogen	antigen sequence, corresponding to amino acids 894-1041 of Human CCDC129.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [CCDC129 coiled-coil domain containing 130 \[Homo sapiens \]](#)

Official Symbol	CCDC129
Synonyms	CCDC129; coiled-coil domain containing 129; coiled-coil domain-containing protein 129;
Entrez Gene ID	223075
Protein Refseq	NP_001244896.1
UniProt ID	A0A087WUB1
Function	receptor binding;