



Anti-NCR3LG1 (aa 38-169) polyclonal antibody (DPABH-20113)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Triggers NCR3-dependent natural killer cell activation.
Immunogen	antigen sequence, corresponding to amino acids 38-169 of Human B7-H6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	NCR3LG1 natural killer cell cytotoxicity receptor 3 ligand 2 [Homo sapiens]
-----------	---

Official Symbol	NCR3LG1
Synonyms	NCR3LG1; natural killer cell cytotoxicity receptor 3 ligand 1; B7H6; B7-H6; DKFZp686O24166; natural cytotoxicity triggering receptor 3 ligand 1; B7 homolog 6; putative Ig-like domain-containing protein DKFZp686O24166/DKFZp686I21167;
Entrez Gene ID	374383
Protein Refseq	NP_001189368.1
UniProt ID	Q68D85
Function	structural molecule activity;