



Anti-KIF18B (aa 15-111) polyclonal antibody (DPABH-20117)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play an important role in regulating cell division.
Immunogen	antigen sequence, corresponding to amino acids 15-111 (STLQVVVRVR PPTPRELDSQ RRPVVQVVDE RVLVFNPEEP DGGFPGLKWG GTHDGPKKKG KDLTFVFDRV FGEAATQQDV FQHTTHSVLD SFLQGYN) of Human KIF18B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 0.02Preservative: 0.02% Sodium azide; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	KIF18B kinesin family member 19B [Homo sapiens]
Official Symbol	KIF18B
Synonyms	KIF18B; kinesin family member 18B; kinesin-like protein KIF18B;
Entrez Gene ID	146909
Protein Refseq	NP_001251503.1
UniProt ID	Q6NWX8
Function	ATP binding; kinesin binding; microtubule binding; microtubule motor activity