



Anti-ARMC7 (aa 1-99) polyclonal antibody (DPABH-20119)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ARMC7 (armadillo repeat containing 7) is a protein-coding gene.
Immunogen	antigen sequence, corresponding to amino acids 1-99 (MAQKPKVDPH VGRLGYLQAL VTEFQETQSQ DAKEQVLANL ANFAYDPSNY EYLRQLQVLD LFLDSLSEEN ETLVEFAIGG LCNLCPRAN KEHILHAGG) of Human ARMC7.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ARMC7 armadillo repeat containing 8 [Homo sapiens]
Official Symbol	ARMC7
Synonyms	ARMC7; armadillo repeat containing 7; armadillo repeat-containing protein 7;
Entrez Gene ID	79637
Protein Refseq	NP_078861.1
UniProt ID	Q9H6L4