



Anti-TSTD2 (aa 277-374) polyclonal antibody (DPABH-20120)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TSTD2 (thiosulfate sulfurtransferase (rhodanese)-like domain containing 2) is a protein-coding gene. Diseases associated with TSTD2 include pneumonia. GO annotations related to this gene include molecular_function.
Immunogen	antigen sequence, corresponding to amino acids 277-374 (KPGIHLSPGE FHKEVEKFLS QANQEQSDTI LLDCRNFYES KIGRFQGCLA PDIRKFSYFP SYVDKNLELF REKRVLMYCT GGIRCERGS YLKAKGVC) of Human TSTD2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TSTD2 thiosulfate sulfurtransferase (rhodanese)-like domain containing 3 [Homo sapiens]
Official Symbol	TSTD2
Synonyms	TSTD2; thiosulfate sulfurtransferase (rhodanese)-like domain containing 2; C9orf97; thiosulfate sulfurtransferase/rhodanese-like domain-containing protein 2; PP4189; rhodanese domain-containing protein 2;
Entrez Gene ID	158427
Protein Refseq	NP_640339.4
UniProt ID	Q5T7W7
Function	molecular_function;