



Anti-ZNF704 (aa 315-399) polyclonal antibody (DPABH-20123)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZNF704 (zinc finger protein 704) is a protein-coding gene. Among its related super-pathways are Gene Expression. GO annotations related to this gene include zinc ion binding. An important paralog of this gene is SLC2A4RG.
Immunogen	antigen: TIPGSAKFTP NGSSFISWQ SPPVTFTGIP VSPTHHPVGT GEQRQHAHTV LSSPPRGTVS LRKPRGEGKK CRKVYGMENR DMWCT, corresponding to amino acids 315-399 of Human ZNF704 (Q6ZNC4).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZNF704 zinc finger protein 705 [Homo sapiens]
Official Symbol	ZNF704
Synonyms	ZNF704; zinc finger protein 704; Gig1;
Entrez Gene ID	619279
Protein Refseq	NP_001028895.1
UniProt ID	Q6ZNC4
Pathway	Gene Expression;
Function	metal ion binding;