



Anti-PPP1R42 (aa 142-223) polyclonal antibody (DPABH-20127)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PPP1R42 (protein phosphatase 1, regulatory subunit 42) is a protein-coding gene. GO annotations related to this gene include actin binding and dynein binding.
Immunogen	antigen sequence: HSLAKSLCIL NISNNNIDDI TDLELLENLN QLIAVDNQLL HVKDLEFLN KLMKLWKIDL NGNPVCLKPK YRDRLILVSK SLGT, corresponding to amino acids 142-223 of Human LRRC67.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PPP1R42 protein phosphatase 1, regulatory subunit 43 [Homo sapiens]
Official Symbol	PPP1R42
Synonyms	PPP1R42; protein phosphatase 1, regulatory subunit 42; dtr; TLLR; LRRC67; protein phosphatase 1 regulatory subunit 42; testis leucine-rich repeat; leucine rich repeat containing 67; leucine-rich repeat-containing protein 67;
Entrez Gene ID	286187
Protein Refseq	NP_001013648.1
UniProt ID	Q7Z4L9
Function	actin binding; dynein binding; tubulin binding;