



Rabbit Anti-Human SMAP2 Polyclonal antibody (DPABH-20130)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SMAP2 fusion protein, sequence: QIQCMQEMGNGKANRLYEAYLPETFRRPQIDPAVEGFIRDKYEKKKYMDRSLDINAFRKE KDDKWKRGSEPVPEKKLEPVVFEKVKMPQKKEDPQLPRKSSPKSTAPVMDLLGLDAPVAC SIANSKTSNTLEKDLDLLASVPSPSSSGSRKVVGSMPTAGSAGSVPENLNLFPPEPGSKSE EIGKKQLSKDSILSLYGSQTPQMPTQAMFMAPAQMAYPTAYPSFPGVTPPNSIMGSMMP PVGMAQPGASGMVAPMAMPAGYMGGMQASMMGVPNGMMTTQQAGYMAGMAAMPQTVYGV QPAQQLQWNLTQMTQQMAGMNFYGANGMMNYGQSMGGNGQAANQTLSPQMWK (77-429aa encoded by BC021133)
------------------	---

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF, ELISA
Positive Control	human liver tissue, human colon tissue, mouse small intestine tissue
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Preservative	0.02% Sodium Azide

Storage

Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Keywords

SMAP2; small ArfGAP2; SMAP1L; stromal membrane-associated protein 2; stromal membrane-associated GTPase-activating protein 2;

GENE INFORMATION

Entrez Gene ID

[64744](#)

UniProt ID

[Q8WU79](#)
