



Anti-ZNF770 (aa 529-618) polyclonal antibody (DPABH-20138)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	antigen sequence: SPYASLCQVE FGNFNNLSNH SGNNVNYNAS QQCQAPGVQK YEVSQDQMS GVKAESQDFI PGSTGQPCLP NVLLESEQSN PFCSYSEHQE, corresponding to amino acids 529-618 of Human ZNF770.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZNF770 zinc finger protein 771 [Homo sapiens]
Official Symbol	ZNF770
Synonyms	ZNF770; zinc finger protein 770; PRO1914;
Entrez Gene ID	54989
Protein Refseq	NP_054825.2
UniProt ID	Q6IQ21
Pathway	Gene Expression.
Function	DNA binding; metal ion binding; sequence-specific DNA binding RNA polymerase II transcription factor activity;