



Anti-B3GNTL1 (aa 9-90) polyclonal antibody (DPABH-20139)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Putative glycosyltransferase.
Immunogen	antigen sequence: ASEESQAMQA HVSIIIPVHN AEPWLDECLR SVLQQDFEGT MELSVFNDAS KDKSGAII EK WRVKLEDSGV HVIIGGHDSP SP, corresponding to amino acids 9-90 of Human B3GNTL1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	B3GNTL1 UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase-like 2 [Homo sapiens]
Official Symbol	B3GNTL1
Synonyms	B3GNTL1; UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase-like 1; B3GNT8; UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase-like protein 1; BGnT-8; beta3Gn-T8; beta3GnTL1; beta-1,3-Gn-T8; BGnT-like protein 1; beta3Gn-T-like protein 1; beta-1,3-N-acetylglucosaminyltransferase 8; beta1,3-N-acetylglucosaminyltransferase-like protein 1; UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase 8;
Entrez Gene ID	146712
Protein Refseq	NP_001009905.1
UniProt ID	A0A024R8X6
Pathway	Metabolism of proteins; O-linked glycosylation of mucins;
Function	transferase activity, transferring glycosyl groups;