



Anti-FAM83H (aa 764-857) polyclonal antibody (DPABH-20146)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a major role in the structural organization and calcification of developing enamel.
Immunogen	antigen sequence: VVSQAWREEV AAPGAVGGER RSLESCLLDL RDSFAQQLHQ EAERQPGAAS LTAAQLLDLTL GRSGSDRLPS RFLSAQSHST SPQGLDSPLP LEGS, corresponding to amino acids 764-857 of Human FAM83H (UniProt ID: Q6ZRV2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	FAM83H family with sequence similarity 84, member H [Homo sapiens]
Official Symbol	FAM83H
Synonyms	FAM83H; family with sequence similarity 83, member H; AI3; protein FAM83H; FAM83H variant 1;
Entrez Gene ID	286077
Protein Refseq	NP_940890.3
UniProt ID	Q6ZRV2