



Anti-SLC38A10 (aa 402-491) polyclonal antibody (DPABH-20148)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Putative sodium-dependent amino acid/proton antiporter.
Immunogen	antigen sequence: PEDLAGEAPG GRLGEAEGLM KVEAARLSAQ DPVVAVAEDG REKPKLPKER EELEQAQIKG PVDVPGREDG KEAPEEAQLD RPGQGIAPPV, corresponding to amino acids 402-491 of Human SLC38A10 (UniProt ID: Q9HBR0).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	SLC38A10 solute carrier family 38, member 11 [Homo sapiens]
Official Symbol	SLC38A10
Synonyms	SLC38A10; solute carrier family 38, member 10; PP1744; putative sodium-coupled neutral amino acid transporter 10;
Entrez Gene ID	124565
Protein Refseq	NP_001033073.1
UniProt ID	Q9HBR0