



Anti-SAMD12 (aa 6-101) polyclonal antibody (DPABH-20149)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SAMD12 (sterile alpha motif domain containing 12) is a protein-coding gene. An important paralog of this gene is SAMD10.
Immunogen	antigen sequence: LHCGLNPRGI DHPAHAEGIK LQIEGEGVES QSIKNKNFQK VPDQKGTPKR LQAEAETAKS ATVKLSKPVA LWTQQDVCKW LKKHCPNQYQ IYSESF, corresponding to amino acids 6-101 of Human SAMD12 (UniProt ID:Q8N8I0).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	SAMD12 sterile alpha motif domain containing 13 [Homo sapiens]
Official Symbol	SAMD12
Synonyms	SAMD12; sterile alpha motif domain containing 12; sterile alpha motif domain-containing protein 12; SAM domain-containing protein 12;
Entrez Gene ID	401474
Protein Refseq	NP_001095146.1
UniProt ID	B8ZZB7