



Anti-RPF1 (aa 11-106) polyclonal antibody (DPABH-20151)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be required for ribosome biogenesis.
Immunogen	antigen sequence: SGKKSLKRKA AAELQEAAG AGDGATENG V QPPKAAAFPP GFSISEIKNK QRRHLMFTRW KQQQRKEKLA AKKKLKKERE ALGDKAPPKP VPKTID, corresponding to N terminal amino acids 11-106 of Human RPF1 (UniProt ID: Q9H9Y2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	RPF1 ribosome production factor 2 homolog (S. cerevisiae) [Homo sapiens]
Official Symbol	RPF1
Synonyms	RPF1; ribosome production factor 1 homolog (S. cerevisiae); BXDC5; ribosome production factor 1; 2310066N05Rik; brix domain containing 5; RNA processing factor 1 (RPF1); brix domain-containing protein 5; ribosome biogenesis protein RPF1;
Entrez Gene ID	80135
Protein Refseq	NP_079341.2
UniProt ID	Q9H9Y2
Function	RNA binding; poly(A) RNA binding; rRNA binding;