



Anti-C9ORF171 (aa 81-160) polyclonal antibody (DPABH-20154)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C9orf171 (chromosome 9 open reading frame 171) is a protein-coding gene.
Immunogen	antigen sequence: SSAVQKVIPS LAGHHIKGGP QAELGKPRER SYSLPGINFN YGLYIRGLDG GVPEAIGRWV VFKQQPTCPH ELTRNYIAMN, corresponding to amino acids 81-160 of Human C9orf171.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C9ORF171 chromosome 9 open reading frame 172 [Homo sapiens]
Official Symbol	C9ORF171
Synonyms	C9ORF171; chromosome 9 open reading frame 171; uncharacterized protein C9orf171;
Entrez Gene ID	389799
Protein Refseq	NP_001269886.1
UniProt ID	Q6ZQR2