



Anti-C8ORF58 (aa 24-98) polyclonal antibody (DPABH-20159)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C8orf58 (chromosome 8 open reading frame 58) is a protein-coding gene. An important paralog of this gene is ENSG00000248235.
Immunogen	antigen sequence: CIVPGVTSTY RRPDAAHGC SSWERGDKFR GVGREALFLK LASRDSGVEM AVGDSPLAAL PGLSQDSLDF ESSGS, corresponding to amino acids 24-98 of Human C8orf58.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C8ORF58 chromosome 8 open reading frame 59 [Homo sapiens]
Official Symbol	C8ORF58
Synonyms	C8ORF58; chromosome 8 open reading frame 58; uncharacterized protein C8orf58;
Entrez Gene ID	541565
Protein Refseq	NP_001013864.1
UniProt ID	Q8NAV2