



Anti-XKR6 (aa 495-591) polyclonal antibody (DPABH-20164)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	XKR6 (XK, Kell blood group complex subunit-related family, member 6) is a protein-coding gene. Diseases associated with XKR6 include keratolytic winter erythema , and hydrocele. An important paralog of this gene is XKR8.
Immunogen	antigen sequence: YGVLHPTGPR AKILASSCCA ELLWGIPLPP DVEPMAPEIP GYRGTVPTPT RAVTEQQEDL TADTCLPVFQ VRPMGPPTPL GRPYLPEGPL IKIDMPR, corresponding to amino acids 495-591 of Human XKR6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	XKR6 XK, Kell blood group complex subunit-related family, member 7 [Homo sapiens]
Official Symbol	XKR6
Synonyms	XKR6; XK, Kell blood group complex subunit-related family, member 6; XRG6; C8orf5; C8orf7; C8orf21; XK-related protein 6; Transmembrane protein C8orf5; X Kell blood group precursor-related family, member 6;
Entrez Gene ID	286046
Protein Refseq	NP_775954.2
UniProt ID	Q5GH73