



Anti-ZNF114 (aa 63-136) polyclonal antibody (DPABH-20167)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	antigen sequence: LPKRTFPEAN RVCLTSSSQ HSTLREDWRC PKTEEPHRQG VNNVKPPAVA PEKDESPVSI CEDHEMRNHS KPTC, corresponding to amino acids 63-136 of Human ZNF114.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZNF114 zinc finger protein 115 [Homo sapiens]
Official Symbol	ZNF114
Synonyms	ZNF114; zinc finger protein 114; zinc finger protein 114 (Y18);
Entrez Gene ID	163071
Protein Refseq	NP_705836.1
UniProt ID	Q8NC26
Pathway	Gene Expression.
Function	DNA binding; metal ion binding; molecular_function;