



Anti-C8ORF86 (aa 152-221) polyclonal antibody (DPABH-20171)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C8orf86 (chromosome 8 open reading frame 86) is a protein-coding gene.
Immunogen	antigen sequence: VLEPIKLFPV SSLRSPLCLN CGSCRESIRI SGELIGNAHS PAPRTPPELE TLGWDKQAVL SGAQVILVCA, corresponding to C terminal amino acids 152-221 of Human C8orf86.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C8ORF86 chromosome 8 open reading frame 87 [Homo sapiens]
Official Symbol	C8ORF86
Synonyms	C8ORF86; chromosome 8 open reading frame 86; uncharacterized protein C8orf86;
Entrez Gene ID	389649
Protein Refseq	NP_997295.1
UniProt ID	Q6ZUL3