



Anti-ERICH5 (aa 204-279) polyclonal antibody (DPABH-20172)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen sequence: QPQGTVGKDE QAPLLETISK ENESPEILEG SQFVETAEED QLQATLGKEE QPQLLERIPK ENVTPEVLDR SQLVEK, corresponding to internal sequence amino acids 204-279 of Human C8orf47.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ERICH5 glutamate-rich 6 \[Homo sapiens \]](#)

Official Symbol	ERICH5
Synonyms	ERICH5; glutamate-rich 5; C8orf47; glutamate-rich protein 5;
Entrez Gene ID	203111
Protein Refseq	NP_001164277.1
UniProt ID	Q6P6B1