



Anti-SPRTN (aa 323-403) polyclonal antibody (DPABH-20173)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SPRTN (SprT-like N-terminal domain) is a protein-coding gene. GO annotations related to this gene include ubiquitin binding and K63-linked polyubiquitin binding.
Immunogen	antigen sequence: SHQNVLSNYF PRVSFANQKA FRGVNGSPRI SVTVGNIPKN SVSSSSQRRV SSSKISLRNS SKVTESASVM PSQDVSGSED T, corresponding to internal sequence amino acids 323-403 of Human C1orf124.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SPRTN SprT-like N-terminal domain [Homo sapiens]
Official Symbol	SPRTN
Synonyms	SPRTN; SprT-like N-terminal domain; DVC1; PRO4323; Spartan; C1orf124; DDDL1880; dJ876B10.3; sprT-like domain-containing protein Spartan; DNA damage protein targeting VCP; SprT-like domain at the N terminus; DNA damage-targeting VCP (p97) adaptor; protein with SprT-like domain at the N terminus; zinc finger RAD18 domain-containing protein C1orf124;
Entrez Gene ID	83932
Protein Refseq	NP_001010984.1
UniProt ID	Q9H040
Function	DNA binding; K63-linked polyubiquitin binding; metal ion binding; protein binding