



Anti-CNPY4 (aa 25-97) polyclonal antibody (DPABH-20176)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Plays a role in the regulation of the cell surface expression of TLR4.
Immunogen	antigen sequence: KEEDDDTERL PSKCEVCKLL STELQAELSR TGRSREVLEL GQVLDTGKRK RHPYYSVSET RLEEALENLC ERI, corresponding to internal sequence amino acids 25-97 of Human CNPY4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CNPY4 canopy FGF signaling regulator 5 [Homo sapiens]
Official Symbol	CNPY4
Synonyms	CNPY4; canopy FGF signaling regulator 4; PRAT4B; protein canopy homolog 4; canopy 4 homolog; PRotein Associated with Tlr4;
Entrez Gene ID	245812
Protein Refseq	NP_689968.1
UniProt ID	Q8N129
Function	receptor binding;