



Anti-ERICH5 (aa 286-368) polyclonal antibody (DPABH-20178)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen antigen sequence: FHKTPEGPGN MEQIQPEGIV GSMEHPARNV EAGAYVEMIR NIHTNEEDQR IEGETGEKVE TDMENEKVSE GAETKEEETG EVV, corresponding to C terminal amino acids 286-368 of Human C8orf47.

Isotype IgG

Source/Host Rabbit

Species Reactivity Human

Purification Immunogen affinity purified

Conjugate Unconjugated

Applications ICC/IF, IHC-P

Format Liquid

Size 100 µl

Buffer pH: 7.20; Constituents: 59% PBS, 40% Glycerol

Preservative 0.02% Sodium Azide

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ERICH5 glutamate-rich 6 \[Homo sapiens \]](#)

Official Symbol	ERICH5
Synonyms	ERICH5; glutamate-rich 5; C8orf47; glutamate-rich protein 5;
Entrez Gene ID	203111
Protein Refseq	NP_001164277.1
UniProt ID	Q6P6B1