



Anti-TMEM242 (aa 1-25) polyclonal antibody (DPABH-20182)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM242 (transmembrane protein 242) is a protein-coding gene.
Immunogen	antigen sequence: METAGAATGQPASGLEAPGSTNDRL, corresponding to N terminal amino acids 1-25 of Human C6orf35.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [TMEM242 transmembrane protein 243 \[Homo sapiens \]](#)

Official Symbol	TMEM242
Synonyms	TMEM242; transmembrane protein 242; BM033; C6orf35; UPF0463 transmembrane protein C6orf35;
Entrez Gene ID	729515
Protein Refseq	NP_060922.2
UniProt ID	Q9NWH2