



Anti-ARMC2 (aa 251-332) polyclonal antibody (DPABH-20183)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ARMC2 (armadillo repeat containing 2) is a protein-coding gene. Diseases associated with ARMC2 include cleft lip , and prostate cancer.
Immunogen	antigen sequence: PKADLQEEDA EIEVDEVFWN TRIVPILREL EKEENIETVC AACTQLHHAL EEGNMLGNKF KGRSILLKTL CKLVDVGSDS LS, corresponding to internal sequence amino acids 251-332 of Human ARMC2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ARMC2 armadillo repeat containing 3 [Homo sapiens]
Official Symbol	ARMC2
Synonyms	ARMC2; armadillo repeat containing 2; bA787I22.1; armadillo repeat-containing protein 2;
Entrez Gene ID	84071
Protein Refseq	NP_001273538.1
UniProt ID	Q8NEN0