



# Anti-C7ORF72 (aa 364-438) polyclonal antibody (DPABH-20184)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | C7orf72 (chromosome 7 open reading frame 72) is a protein-coding gene. GO annotations related to this gene include molecular_function.                         |
| <b>Immunogen</b>           | antigen sequence: ANIPGYTGKV HFTATHPANS NIPSTTPSPD SELHRVFQKE<br>MAVDLFRHQA PLSRLVTTVR PYNPFNKKDK ETIDY, corresponding to amino acids 364-438 of Human C7orf72 |
| <b>Isotype</b>             | IgG                                                                                                                                                            |
| <b>Source/Host</b>         | Rabbit                                                                                                                                                         |
| <b>Species Reactivity</b>  | Human                                                                                                                                                          |
| <b>Purification</b>        | Immunogen affinity purified                                                                                                                                    |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                   |
| <b>Applications</b>        | ICC/IF, IHC-P                                                                                                                                                  |
| <b>Format</b>              | Liquid                                                                                                                                                         |
| <b>Size</b>                | 100 µl                                                                                                                                                         |
| <b>Buffer</b>              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                                                                                                                  |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                                                                             |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.                                                                          |

## GENE INFORMATION

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">C7ORF72 chromosome 7 open reading frame 73 [ Homo sapiens ]</a>   |
| <b>Official Symbol</b> | C7ORF72                                                                       |
| <b>Synonyms</b>        | C7ORF72; chromosome 7 open reading frame 72; uncharacterized protein C7orf72; |
| <b>Entrez Gene ID</b>  | <a href="#">100130988</a>                                                     |
| <b>Protein Refseq</b>  | <a href="#">NP_001155306.2</a>                                                |
| <b>Function</b>        | molecular_function;                                                           |