



Rabbit Anti-Human DENND3 Polyclonal antibody (DPABH-20185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DENND3 fusion protein, sequence: MRSLRKKREKPRPEQWKGLPGPPRAPEPEDVAVPGGVDLLTLPQLCFPGGVCVATEPKED CVHFLVLTDCGNRTYGVVAQYYRPLHDEYCFYNGKTHRECPGCFVPFAVCVVSFRFPYYN SLKDCLSCLLALLKPKCDFEVD SHIKDFAAKLSLIPSPPPGPLHLVFNMKSLQIVLPARA DPESPILDLDLHPLLCFRPEKVLQILTCILTEQRIVFFSSDWALLTLVTECFMAYLYPL QWQHFPFVPI LSDQMLDFVMAPTSFLMGCHLDHFEEVSKEADGLVLINIDHGSITYSKSTD DNVDIPDVPLLA AQTFIQRDVKNHNLN YEHRVFNSEEF LKTRAPGDHQFYKQVLD TYMFHS (1-360 aa encoded by BC133023)
------------------	---

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Positive Control	mouse stomach tissue
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Preservative	0.1% Sodium Azide

Storage

Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Introduction

DENND3 may play a role in protein transport from recycling endosomes to lysosomes regulating, for instance, the degradation of the transferrin receptor

Keywords

DENND3; DENN/MADD domain containing 3; DENN domain-containing protein 3;

GENE INFORMATION

Entrez Gene ID

[22898](#)

UniProt ID

[A2RUS2](#)
