



Rabbit anti-Human ARHGEF10L polyclonal antibody (DPABH-20187)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Acts as guanine nucleotide exchange factor (GEF) for RHOA, RHOB and RHOC.
Immunogen	antigen sequence: ALEHSEEDGS IYEMADDPDV WVRSRPCARD AHRKEICSV IISGGQGYRN FGSALGSSGR QAPCGETDST LLIWQVPL, corresponding to C terminal amino acids 1200-1277 of Human ARHGEF10L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ARHGEF10L Rho guanine nucleotide exchange factor (GEF) 11-like [Homo sapiens]
Official Symbol	ARHGEF10L
Synonyms	ARHGEF10L; Rho guanine nucleotide exchange factor (GEF) 10-like; GrinchGEF; rho guanine nucleotide exchange factor 10-like protein;
Entrez Gene ID	55160
Protein Refseq	NP_001011722.2
UniProt ID	Q9HCE6
Pathway	Regulation of RhoA activity;
Function	Rho guanyl-nucleotide exchange factor activity;